



Labflex laboratory chairs

Sophisticated seating solutions for your laboratory

New generation of lab chairs

Laboratory chair for all kinds of laboratories. The lab chair offers ergonomics with adaptability, comfort and intuition.

WORK IS CHANGING

Today's complex working environment becomes more and more demanding. A company that wants to attract and keep qualified employees has to provide ideal working conditions. These ideal working conditions start with an ideal workplace chair.



COMFORTABLE

The ergonomic upholstery offers the best possible support for the body, outstanding sitting climate and top comfort.



ADJUSTMENT

The intelligently designed operational controls facilitate perfect adjustment for any user.


Laboratory chairs

New generation of workplace chairs



73349573	73349571
Seat height 460 – 620 mm	Seat height 590 – 870 mm
With casters	With glides and footrest
Backrest height 570 mm	Backrest height 570 mm

- Ergonomic and flexible with weight adjustment, seat depth adjustment, adjustable seat inclination and backrest adjustable in height.
- Excellent seating comfort.
- Washable and resistant to disinfectants.
- Easy to clean thanks to surface with as few joints as possible.
- Sustainability and flexibility thanks to changeable upholstery elements.



Fraunhofer
IAO

Basic laboratory chairs

Laboratory chairs for all kinds of laboratories. The perfect combination of comfort and ergonomics.




3D-FLEX FUNCTION OF THE BACKREST

The backrest adapts to the movements of the user following the upper body body when leaning back or turning sideways.



COMFORT AND ERGONOMICS

The front edge of the seat automatically follows the forward-angled working position of the body. This ensures optimal blood circulation and effectively prevents pressure pain and dead legs.



Basic laboratory chairs

Laboratory chairs for all kinds of laboratories. The perfect combination of comfort and ergonomics.





73349123	73349125	73349127
Seat height 450 – 650 mm	Seat height 560 – 810 mm	Seat height 450 – 650 mm
With casters	With Stop & Go casters	With casters
Backrest height 420 mm	Backrest height 420 mm	Large, comfortable seat

- Hygienic design with minimal joints for easy cleaning and disinfection.
- Resistant to all usual disinfectants and chemicals in accordance with ISO 2812.
- Antibacterial upholstered covering.
- Suitable for cleanrooms category 3 in accordance with ISO 14644-1.
- Adheres to biosafety level BSL1, BSL2 and BSL3 of the Directive of Biological Agents.



Basic laboratory chair

The tried-and-tested all-rounder for use in laboratories.



73349133

Seat height 470 - 610 mm

With casters

Backrest height 420 mm

Backrest adjustable in height

Steel base

Resistant to disinfectants



Garantie
10
Jahre

Cleanroom laboratory work

Laboratory chairs suitable for use in cleanrooms and designed to meet the specifications of the GMP regulations.



73349103	73349101	73349107
Seat height 400 – 510 mm	Seat height 550 – 800 mm	Seat height 450 – 650 mm
With casters	With glides	With casters
Backrest height 420 mm	Backrest height 420 mm	Large, comfortable seat

- Auto-Motion technology supports the user's posture
- Suitable for cleanrooms in air purity class 3 in accordance with DIN EN ISO 14644-1.
- Meets the specifications of the GMP regulations.
- No room for dust and dirt thanks to seamless design.
- Washable and resistant to disinfectants.
- Mechanism with integrated control panel protected by soft cover.
- Compact and space saving.



Laboratory standing rests

Laboratory standing rests are the best support for sitting, standing and combined workplaces. Our standing rests meet the ergonomic requirements of standing workstations.



73349144	73349408
Seat height 620 – 850 mm	Seat height 450 – 650 mm
Changeable glides for hard/soft floors	With casters
Height adjustable seat	Height adjustable seat
Comfort seat in Soft Touch PU foam. The seat can be revolved by 20° to the left and to the right	Ergonomic, soft and comfortable seat with back support
Integrated seat handles make it easy to handle and store	Seat is 360° turnable
Washable and resistant to disinfectants	Washable and resistant to disinfectants
Easy to clean thanks to surface with as few joints as possible	Easy to clean
Sustainable and flexible use	Usable as both chair and as standing rest



Want to know more or try one of our
lab chairs at your own laboratory?

Contact your local Labflex office.

Labflex is one of Europe's leading suppliers of innovative customized laboratory solutions to industrial and university laboratories and to the research and hospital sector. Labflex offers consulting, project planning, concept development, project management and a wide, well-documented product range as well as installation and service. Further information on www.labflex.com

HEAD OFFICE, COPENHAGEN

Labflex A/S
Marielundvej 45 A
DK - 2730 Herlev

+45 8747 2700

DENMARK

Labflex A/S
Hjortevej 3
DK - 7800 Skive

+45 8747 2700

UNITED KINGDOM

Labflex Limited
Pegasus Business Park
Castle Donington
UK - DE74 2TZ

+44 (0) 1332 638 071